

CHARCOAL ALUMINA COLUMN

● Intended use

The Puritox^{SR} TC-C210 Charcoal Alumina Columns Can Be Used As A Cleanup For Trichothecenes, Including DON Prior To GC/HPLC/TLC Methodology. This Column Can Also Be Used In Conjunction With Other Purification Columns To Enhance The Cleanup Of Complex Matrices.



● General

PuriTox^{SR} Columns are formulated and designed to provide optimum sample extract purification prior to the analyses of mycotoxins by HPLC, GC & TLC methodology. PuriTox^{SR} columns utilize an easy to use one-step purification format that allows you to quickly purify the most complex matrices.

Specifications:

Format:	- 50 Purification Columns
Product ID #:	- TC-C210