

Extracton Supplies - Trilogy Analytical Laboratory - Puritox^{SR} Columns

FUMONISIN COLUMN

● Intended use

The Puritox^{SR} TC-F120 Fumonisin Columns Are Utilized For The Clean Up Of Grains Prior To HPLC Analysis. The TC-F120 Is A Gravity Flow Column That May Be Used With Vacuum.



● General

PuriTox^{SR} Columns are formulated and designed to provide optimum sample extract purification prior to the analyses of mycotoxins by HPLC, GC & TLC methodology. PuriTox^{SR} columns utilize an easy to use one-step purification format that allows you to quickly purify the most complex matrices. PuriTox^{SR} have several advantages over other purification columns.

Specifications:

Format:	- 50 Purification Columns.
Product ID #:	- TC-F120.