

Extraction Supplies - Trilogy Analytical Laboratory - Puritox^{SR} Columns

OCHRATOXIN COLUMN

● Intended use

The Puritox^{SR} TC-O130 Ochratoxin Columns Are Designed For Cleanup Of Grain And Feed Samples Prior To HPLC Or TLC Procedures. These Columns Are Utilized With A Flow Through System With The Aid Of A Vacuum Enhanced Gravity Drain. The Large Reservoir Allows For The Addition Of Increased Volumes Of Extract To Be Placed On The Column.



● General

PuriTox^{SR} Columns are formulated and designed to provide optimum sample extract purification prior to the analyses of mycotoxins by HPLC, GC & TLC methodology. PuriTox^{SR} columns utilize an easy to use one-step purification format that allows you to quickly purify the most complex matrices.

Specifications:

Format:	- 50 purification columns.
Product ID #:	- TC-O130.