

## Extraction Supplies - Trilogy Analytical Laboratory - Puritox<sup>SR</sup> Columns

### PRELISA COLUMN

#### ● Intended use

The Puritox<sup>SR</sup> TC-P500 PreLISA Column Is Ideal For Use On Complex Matrices Prior To ELISA Test Kits. The PreLISA Column Can Be Used For Aflatoxin, DON And Zearalenone To Enhance The Performance Of Elisa Test Procedures.

#### ● General

PuriTox<sup>SR</sup> Columns are formulated and designed to provide optimum sample extract purification prior to the analyses of mycotoxins by HPLC, GC & TLC methodology. PuriTox<sup>SR</sup> columns utilize an easy to use one-step purification format that allows you to quickly purify the most complex matrices.



#### Specifications:

Format:	- 50 Purification Columns
Product ID #:	- TC-P500