

Extraction Supplies - Trilogy Analytical Laboratory - Puritox^{SR} Columns

DON (VOMITOXIN) COLUMN

● Intended use

The PuriTox^{SR} DON Columns TC-T200 are specifically designed for the purification of DON (Vomitoxin) prior to HPLC / GC and other detection methods. These columns will purify grain and feed Sample extracts in one step.

● General

PuriTox^{SR} Columns are formulated and designed to provide optimum sample extract purification prior to the analyses of mycotoxins by HPLC, GC & TLC methodology. PuriTox^{SR} columns utilize an easy to use one-step purification format that allows you to quickly purify the most complex matrices. PuriTox^{SR} have several advantages over other purification columns.



Specifications:

Format:	- 50 Purification Columns.
Product ID #:	- TC-T200.